



Technical Data Sheet

Product Reference: AT19 Autoclave Masking Tape

Description

A masking tape based on smooth crepe paper and rubber adhesive

The product is designed as an indicator for steam autoclaving (steam sterilisation)

During service black lines will appear to indicate when required steam sterilisation temperature has been attained

Properties

Carrier	:	Creped Paper
Adhesive	:	Rubber
Thickness	:	0.17mm
Tensile	:	90N/25mm
Adhesion to Steel	:	10N/25mm
Elongation	:	7%
Temperature Range	:	125°C to 135°C / 180 sec to 240 sec
Shelf Life stored	:	1 year in original packaging
Product complies with Specifications	:	in cool dry conditions UU-T-118B (US) C.I.D. A-A 54589A (US)

Statements, technical information and recommendations contained herein are based on production data from the manufacturer with properties quoted being typical but intended only as a guide. They do not constitute a product specification nor are they to be construed in any manner as warranties expressed or implied. The User shall determine the suitability of the product for its intended use and the User assumed all risk and liability whatsoever in connection therewith.

The Old Bakery, Ashley Gardens, Tunbridge Wells, Kent,
TN4 8TU, United Kingdom

Tel: +44 (0)1892 514444 • **USA call tollfree** to 1 877 412 4876

Fax: +44 (0)1892 514411

Email: sales@temperature-indicators.co.uk