

## Autoclave Tape ISO9000 approved

### Properties

Carrier:	Creped
paper Adhesive:	Rubber
Thickness:	0.185mm
Tensile:	93N/25mm
Adhesion to steel:	11N/25mm
Elongation:	7%
Temperature range:	125°C - 135°C
	180sec to 240sec
Shelf life:	1 year (original
packaging) Storage:	Cool, dry conditions

### Description

- A tape based on a smooth crepe paper and rubber adhesive.
- The product is designed as an indicator for steam autoclaving. (Steam sterilisation).
- During service black lines will appear to indicate when steam sterilisation has been attained.

### Supplied

- 19mm wide - 48 rolls per box
- 25mm wide - 24 rolls per box
- 50mm wide - 24 rolls per box



Please telephone **01892 514444**  
for your **FREE** sample pack